



BS-3014 Stereo Microscope

Instruction Manual



BS-3014A



BS-3014B



BS-3014C



BS-3014D

| | |
|---|-----------|
| 1. Before use | 1 |
| 2. Nomenclature | 2 |
| 3. Assemblage | 3 |
| 4. Operation | |
| 4-1 Use the glass stage..... | 6 |
| 4-2 Adjust the degree of tightness of the focusing arm..... | 6 |
| 4-3 Set the specimen side..... | 6 |
| 4-4 Adjust diopter and focus..... | 6 |
| 4-5 Adjust the interpupillary distance..... | 6 |
| 4-6 Use Eyepiece shields..... | 7 |
| 4-8 Install the illumination device..... | 7 |
| 4-9 Choose the optical system..... | 7 |
| 4-10 Mount the photo eyepiece and the PK-mount adapter..... | 7 |
| 4-11 Adjust the brightness of the bottom light..... | 8 |
| 4-12 Replace the lamps..... | 8 |
| 4-13 Replace the fuse..... | 8 |
| 5. Configuration | 9 |
| 6. Technical parameter | 10 |
| 7. Troubleshooting | 12 |

1 Before use


1-1 **Notice**


- 1) Microscope ought to be placed in a dry and clean place. Do not expose the microscope in the sun directly. Avoid high temperature and violent vibration.
- 2) Microscope is a precision instrument, so handle with care, avoiding impact or abrupt movement during transportation.
- 3) To keep the image clear, do not leave fingerprints or stains on the surfaces of the lens.
- 4) Never turn the left and right focusing knob in the adverse direction at the same time, otherwise the microscope will be damaged.

1-2 **Maintenance**


- 1) All lenses must be kept clean. Fine dust on the surface of the lens should be blown off with a hand blower or wiped off gently with a soft lens tissue; Fingerprints or oil marked on it should be wiped off with a tissue moistened with a small amount of a 3:7 mixture of alcohol and ether.
- 2) Never use the organic solution to clean the other surface (especially the plastic surfaces). If necessary, please choose the neutral detergent.
- 3) Do not take the microscope apart for fearing that it is damaged.
- 4) After using, cover the microscope with the dust-cover provided and store it in a dry and clean place free from moisture to prevent rust.
- 5) To keep the performance of the microscope, please check it periodically.


2 Nomenclature





3 Assemblage





4 Operation




Fig. 1




Fig. 2





Fig. 3



4-1 Use the glass stage

- 1) Press the glass stage on the sunken place then the other side of the glass stage will be lifted. (Fig. 1)

4-2 Adjust the degree of tightness of the focusing arm

- 1) If you want to adjust the degree of tightness of the focusing arm, you can hold one of the focusing knob and turn another one to attain a suitable position. The degree of tightness relies on the direction to be turned. The clockwise direction is tight, otherwise, is loose.
- 2) The suitable position of the tightness can make the adjustment more comfortable and prevent the focusing bracket from slipping down by its weight during the observation. (Fig. 2)

4-3 Set the specimen slide

- 1) Set the specimen in the center of stage plate. If necessary, clamp the slide with the clips.
- 2) Turn on the light.

4-4 Adjust diopter and focus

- 1) Turn the focusing knob and observe the specimen through the right eyepiece till the image of the specimen is clear.
- 2) Observe the specimen through the left eyepiece and adjust the diopter adjustment ring ① till the image is clear. (Fig. 3)

4-5 Adjust the interpupillary distance

- Adjust the prism housing along the direction of arrow of the Fig. 4 till the observation is comfortable.

4-6 Use eyepiece shields

- 1) For user who does not wear glasses, hold the diopter-adjusting ring to prevent them from rotating and turn the eyepiece till the eyepiece shield fit the observer well.
- 2) For user who wears glasses, take the eyepiece shields off before observation.




Fig. 5

4-7 Install and remove the optional eyepiece micrometer

- 1) Turn and remove the mounting ring (2) from the eyepiece. (Fig. 5)
- 2) Clean the eyepiece micrometer (1), and mount it to the mounting ring with the inscription side downward.
- 3) Gently twist the mounting ring with the eyepiece micrometer into the eyepiece till tightening (2) securely.
- 4) To remove the eyepiece micrometer, take down the mounting ring by twisting and take out of the micrometer, and wrap it in clean soft paper for storage.




Fig. 6

4-8 Install the illumination device

- 1) Insert the illumination device (1) in the bracket with the protrudent side toward the lock-screw (2) and tighten the lock-screw. (Fig. 6)
- 2) Put the plug into the socket of the pillar stand (3).




Fig. 7

4-9 Choose the optical system

- 1) You can alternate the binocular and video capture by pushing or pulling "the pole". You can attain binocular observation by pushing "the pole" inside, or attain video capture by pulling it outside. (Fig. 7)

4-10 Mount the photo eyepiece and the PK-mount adapter




Fig. 8




Fig. 9

- 1) Put the photo eyepiece into the eyepiece socket of the tri-ocular.
- 2) Connect the PK-mount adapter with the photo eyepiece, and then tighten the lock-screw. (Fig.8)

4-11 Adjust the brightness of the bottom light

- 1) Turn the adjustable light knob ① according to the sign marked on the base, along the clockwise the brightness will be added, otherwise it will be weakened. (Fig.9)




Fig.10

4-12 Replace the lamps

- 1) Press the stage on the sunken place then the other side will be lifted. (Fig.10)
- 2) Take the lamp out of the jack.
- 3) Put a new lamp into the jack thoroughly.
- 4) Recover the stage plate. (Fig.11)

Note: ① Before replacing the lamps, turn off the power first.

② Avoid violence while the lamp is plugged into the jack.




Fig.11

4-13 Replace the fuse

- 1) Screw the fuse tube out with a screwdriver and then pull the fuse out of the tube ①.
- 2) Remove the fuse and mount it in an adverse way. (Fig. 12)




Fig. 12

5 Configuration chart

5-1 BS-3014 Series Configuration

| Item | Specification | BS-3014A | BS-3014B | BS-3014C | BS-3014D |
|---------------------|---|----------|----------|----------|----------|
| Head | Binocular Viewing Head, Inclined at 45°, 360° rotatable, Interpupillary adjusting distance 54-76mm, left eyepiece with diopter adjustment ± 5 | ● | ● | ● | ● |
| Eyepiece | High eyepoint WF10 \times /20mm eyepiece | ● | ● | ● | ● |
| | WF15 \times /15mm eyepiece | ○ | ○ | ○ | ○ |
| | WF20 \times /10mm eyepiece | ○ | ○ | ○ | ○ |
| Objective | 2 \times , 4 \times | ● | ● | ● | ● |
| | 1 \times , 2 \times | ○ | ○ | ○ | ○ |
| | 1 \times , 3 \times | ○ | ○ | ○ | ○ |
| Magnification | 20 \times , 40 \times , with optional eyepiece and auxiliary objective, can be extended to 5 \times -160 \times | ● | ● | ● | ● |
| Auxiliary Objective | 0.5 \times objective, W.D.: 165mm | ○ | ○ | ○ | ○ |
| | 1.5 \times objective, W.D.: 45mm | ○ | ○ | ○ | ○ |
| | 2 \times objective, W.D.: 30mm | ○ | ○ | ○ | ○ |
| Working Distance | 100mm | ● | ● | ● | ● |
| Head Mount | 76mm | ● | ● | ● | ● |
| Illumination | Transmitted light 12V/15W Halogen, Brightness Adjustable | | ● | | |
| | Incident light 12V/15W Halogen, Brightness Adjustable | | ● | | |
| | Transmitted light 3W LED, Brightness Adjustable | | ○ | | ● |
| | Incident light 3W LED, Brightness Adjustable | | ○ | | ● |
| | LED ring light | ○ | ○ | ○ | ○ |

| | | | | | |
|--------------|---|---|---|---|---|
| | Cold light source | ○ | ○ | ○ | ○ |
| Focusing Arm | Coarse focusing, focusing range 50mm | ● | ● | ● | ● |
| Pillar Stand | Pole height 240mm, pole diameter Φ32mm, with Clips, Φ95 black&White plate, Base size: 200×255×22mm, no illumination | ● | | | |
| | Pole height 240mm, pole diameter Φ32mm, with Clips, Φ95 black&White plate, glass plate, Base size: 200×255×60mm, Halogen illumination | | ● | | |
| | Pole height 240mm, pole diameter Φ32mm, with Clips, Φ95 black&White plate, Base size: 205×275×22mm, no illumination | | | ● | |
| | Pole height 240mm, pole diameter Φ32mm, with Clips, Φ95 black&White plate, glass plate, Base size: 205×275×40mm, LED illumination | | | | ● |
| Package | 1pc/1carton, 38.5cm*24cm*37cm, Net/Gross Weight: 3.5/4.5kg | ● | ● | ● | ● |

Note: ● Standard Outfit, ○ Optional

6 Technical parameter

6-1 BS-3014 series optical parameter

| Objective Mag. | Working Distance (mm) | Eyepiece | | Eyepieces (option) | | | |
|-------------------|-----------------------------|------------|-----------------|--------------------|-----------------|------------|-----------------|
| | | WF10X/20mm | | WF15X/15mm | | WF20X/10mm | |
| | | Mag. | Objective field | Mag. | Objective field | Mag. | Objective field |
| 1X | 100 | 10X | 20 | 15X | 15 | 20X | 11 |
| 2X | | 20X | 10 | 30X | 7.5 | 40X | 5 |
| 3X | | 30X | 6.7 | 45X | 5 | 60X | 3.3 |
| 4X | | 40X | 5 | 60X | 3.75 | 80X | 2.5 |

6-2 Auxiliary objective for BS-3014 series

| Auxiliary objectives | Magnification | Working distance (mm) |
|----------------------|---------------|-----------------------|
| 0.5 X | 0.5X | 165 |
| 1.5 X | 1.5X | 45 |
| 2 X | 2X | 30 |

★ Working distance is fixed regardless of the magnification factor.

★ Total mag.=Objective mag. X Auxiliary mag. X Auxiliary mag.

Eyepiece field

Diameter of field of view (mm) = $\frac{\text{Eyepiece field}}{\text{Objective mag. X Auxiliary objective mag.}}$

★ Photo adaptor mag.=Objective mag.(X Auxiliary objective mag.)X Photo eyepiece mag.

6-3 The base electrical specification of BS-3014 series

| Model | | BS-3014A/C | BS-3014B | BS-3014D |
|--------------|--------------|------------|--------------------------------------|--------------------------------------|
| Parts | | | | |
| Power supply | | No | 220V-50Hz 110V-50/60Hz | 220V-50Hz 110V-50/60Hz |
| Transformer | | No | Input:220/110VAC Output:12VDC/45W | Input:220/110VAC Output:12VDC/45W |
| Illuminate | Top light | No | 12V/15W halogen lamp | 3W LED lamp |
| | Bottom light | | 12V/15W halogen lamp | 3W LED lamp |

6-4

Configuration parameter of BS-3014 series

| Parts \ Model | | BS-3014(1X,2X) | BS-3014(1X,3X) | BS-3014(2X,4X) |
|---------------|------------------------------------|--|----------------|----------------|
| Head | Objective magnification | 1X, 2X | 1X, 3X | 2X, 4X |
| | Working distance | 100mm | | |
| | Observation angle | 45° | | |
| | Interpupillary distance adjustment | Linkage between left and right eyepiece tube Range of single adjustment: 54-75mm | | |
| | Diopter adjustment | Range of single adjustment : ±5D | | |
| | Mount with auxiliary objectives | Screw hole : M48*0.75 | | |
| Objective | Field of view | φ20mm | | |
| Main body | Mount the head | Mount the head in the bracket hole whose diameter is φ76mm | | |
| | Focusing device | The degree of adjustable by rotating the focusing knob. Range of single adjustable :49 mm | | |
| | Glass stage | Diameter: φ95mm | | |
| | Clips | Put it on the base from top | | |

7 Trouble shooting

The performance of the microscope can't be made fully because of unfamiliar using. This table will give some advices.

| Trouble | Cause | Remedy |
|---|---|--|
| 1.Double images | Interpupillary distance is not correct | Readjust it |
| | Diopter adjustment is not correct | Readjust it |
| | Magnification of each eyepiece is not the same size | Mount the same size eyepiece |
| 2.Dirt appears in the field of view | Dirt on the specimen | Clean the specimen |
| | Dirt on the surface of eyepiece | Clean the surface |
| 3.Image is not clear | Dirt on the surface of the objective | Clean the objectives |
| 4.Image is not clear while the focus changing | Diopter adjustment is not correct | Readjust the diopter |
| | Focus is not correct | Readjust the focus |
| 5.The focusing knob is not smooth | The focusing knob is too tight | Loosen it to a suitable position |
| 6.The image is obscure because of the head slipping down by itself during observation | The focusing knob is too loose | Tighten it to a suitable position |
| 7.Incision image appears in the field of view or of the video view | The pole is not in correct position | Pull or push it to the correct position |
| 8.Eyes fell tired easily | Diopter adjustment is not correct | Adjust the diopter |
| | Brightness of light is not correct | Adjust the brightness |
| 9.Bulb does not work when the switch is on | No power in | Check the connection with the power supply |
| | The bulb was not insert correct | Insert it correctly |
| | Bulb is wrong | Replace with a new one |
| 10.Bulb is burned out suddenly | Use the wrong bulb | Replace with a correct one |
| | The voltage is too high | Control the voltage Eg :use voltage regulator |
| 11.Brightness is not enough | Use the wrong bulb | Replace with a correct one |
| | The voltage is too low | Increase the input voltage |
| 12.The bulb flickers or the brightness is unstable | The bulb will burn out soon | Replace with a new one |
| | The bulb was not inserted correctly | Insert it correctly |