

LED-56A / 64A / 144A LED Ring Light



LED-56A



LED-56A(Rotatable)



LED-64A



LED-144A

Application

LED-56A/56A(rotatable)/64A/144A LED ring light are simple and compact, they can be used as incident illumination on stereo microscopes and machine vision systems. LED-56A has a rotatable ring, so easy to change the direction of the LED lamp, these LED ring lights have long working life and save energy features.

Feature

1. Light intensity can be adjusted to get best contrast.
2. LED-56A(rotatable) has a rotatable ring, easy to change the direction of the LED lamp.
3. Reliability and Durability.
4. Low Power consumption, Super Brightness.
5. Reasonable Design, Uniformity Light, Non-flicker & Shadow.

Specification

Model	LED-56A	LED-56A(Rotatable)	LED-64A	LED-144A
LED Quantity	56pcs LED lamps	56pcs LED lamps	64 LED lamps	144 LED lamps
Color	White	White	White	White
Mounting Diameter	Φ61mm	Φ58mm/Φ63mm	Φ61mm	Φ61mm
Light Head Material	Plastic	Plastic	Aluminum	Plastic
Input Power	DC12V, 2W	DC12V, 2W	DC12V, 2.2W	DC12V, 4.2W
Input Voltage	Universal 100-240V AC	Universal 100-240V AC	Universal 100-240V AC	Universal 100-240V AC
Color Temperature	8000K	8000K	8000K	8000K
LED Luminance	0-12000 Lux	0-12000 Lux	0-12000 Lux	0-14000 Lux
LED Lifetime	50,000 hours	50,000 hours	50,000 hours	50,000 hours
Brightness	0-100% adjustable	0-100% adjustable	0-100% adjustable	0-100% adjustable
Working Distance	30-150mm	30-150mm	30-150mm	30-150mm
Standard Set	LED-56A LED ring light, DC12V Power Supply	LED-56A LED ring light, DC12V Power Supply	LED-64A LED ring light, DC12V Power Supply	LED-144A LED ring light, DC12V Power Supply