

## BAL-60A / 60B / 64 LED Ring Light



BAL-60A



BAL-60B/72B



BAL-64

### Application

BAL-60A/60B/72B/64 LED ring light are simple and compact, they can be used as incident illumination on stereo microscopes and machine vision systems. BAL-60B has four-zone lighting control separately, these LED ring lights have long working life and save energy features.

### Feature

1. Light intensity can be adjusted to get best contrast.
2. BAL-60B has Four-zone lighting control separately to enhance contrast and reduce reflection.
3. Reliability and Durability.
4. Low Power consumption, Super Brightness.
5. Reasonable Design, Uniformity Light, Non-flicker & Shadow.

### Specification

Model	BAL-60A	BAL-60B/72B	BAL-64
LED Quantity	60 LED lamps	60/72 LED lamps	64 LED lamps
Color	Black	Black	White
Mounting Diameter	Φ60mm	Φ60mm	Φ60mm
Light Head Material	Plastic	Aluminum	Plastic
Input Power	DC12V	DC12V	DC12V
Input Voltage	Universal 100-240V AC	Universal 100-240V AC	Universal 100-240V AC
Color Temperature	6500K-7000K	6000K-6500K	6500K-7000K
LED Luminance	0-12000 Lux	0-16000 Lux	0-12000 Lux
LED Lifetime	50,000 hours	50,000 hours	50,000 hours
Brightness	0-100% adjustable	0-100% adjustable, 4 zone control separately	0-100% adjustable
Working Distance	30-120mm	60-120mm	30-120mm
Standard Set	BAL-60A LED ring light, DC12V Power Supply	BAL-60B/72B LED ring light, DC12V Power Supply	BAL-60B LED ring light, DC12V Power Supply