

BAL-2, BAL-3 Fluorescent Ring Light



BAL-2B



BAL-2A/2C



BAL-3A



BAL-3B



BAL-3C

With high brightness and even illumination, simple structure and easy to operate, BAL-2, BAL-3 Series fluorescent ring light can be used as incident illumination for various stereo microscopes. The difference of BAL-2A and 2C is the lamp is different.

Model	BAL-2A	BAL-2B	BAL-2C	BAL-3A	BAL-3B	BAL-3C
Voltage	110V/220V					
Power	8w			10w		
Color Temperature	6400K-7000K					
Mounting Diameter	Φ30-Φ60mm					
Light adjustment	No Light adjustment					Light adjustment