

BAL-72 LED Ring Light



Application

BAL-72 LED ring light can be used as incident illumination on stereo microscopes and machine vision systems. It has four-zone lighting control separately, long working life and save energy features.

Feature

1. Intensity Adjustment, Good Variable Control Permit Gradual Increase in Brightness with Four-zone Lighting Control Separately.
2. Reliability and Durability.
3. Low Power consumption, Super Brightness.
4. Metal Shell and ESD.
5. Reasonable Design, Uniformity Light, Non-flicker & Shadow.

Specification

LED Quantity	72LED lamps
Mounting Diameter	Φ62mm
Input Power	DC12V 600mA
Input Voltage	Universal 100-240V AC
Output Power	7.2W
Color Temperature	6500K
LED Luminance	> 15000mcd
LED Lifetime	50,000 hours
Working Distance	95-115mm
Standard Set	BAL-72 LED ring light, Light Control Box, DC12V Power Supply