


BAL2B Series LED Ring Light



BAL2B-60/78

BAL2B series LED ring light have features of high brightness, low temperature and flash free, they can be used as auxiliary illumination for industrial monocular microscopes, stereo microscopes and similar lens.

Feature

1. The light head adopt ABS plastic material and $\phi 5\text{mm}$ LED lamps, with excellent light focus effect and high efficiency.
2. The power controlling adapter adopt aluminum case, solid and reliable.
3. Continues light intensity adjustment can meet different requirement.
4. Reliable circuit board ensures safety and long working life.
5. ESD light head treatment is optional.

Specification

Model	BAL2B-60	BAL2B-78
Input Voltage	Universal 100-240V AC	Universal 100-240V AC
Input Power	5 W	6 W
Mounting Diameter	$\phi 61\text{mm}$	$\phi 70\text{mm}$
LED Quantity	60pcs LED lamps	78pcs LED lamps
LED Lifetime	50,000hrs	50,000hrs
LED Color	White (Other Colors can be customized)	White (Other Colors can be customized)
Color Temperature	6400K, other color temperature can be customized	6400K, other color temperature can be customized
Illuminance@100mm	24000lx	25000lx
Light Control	Brightness Adjustable	Brightness Adjustable
Light head material	ABS Plastic	ABS Plastic
Packing	BAL2B-60 LED ring light head, Aluminum Light Control Box, Power cable	BAL2B-78 LED ring light head, Aluminum Light Control Box, Power cable