

BS-2000A, B, C Biological Microscope



BS-2000A

BS-2000B

BS-2000C

Introduction

With sharp image, competitive and reasonable unit price, BS-2000A, B, C series microscopes are ideal instruments for student use. These microscopes are mainly used in elementary schools.

Specification

Item	Specification	BS-2000A	BS-2000B	BS-2000C
Eyepiece	Huygens Eyepiece H5×, H10×	●	●	●
	Huygens Eyepiece H12.5×	●		
	Huygens Eyepiece H16×		●	●
Objective	Achromatic Objective 10×, 40×	●	●	●
	Achromatic Objective 4×, 60×	○	○	○
	Achromatic Objective 100×	○	○	●
Nosepiece	Triple Nosepiece	●	●	●
Stage	Plain Stage with Slide Clips 110×120mm	●	●	
	Attachable Mechanical Stage 110×120mm/ 60×30mm	○	○	●
Condenser	Single Lens with Disc Diaphragm	●		
	Single Lens with Iris Diaphragm		●	
	Abbe Condenser with Iris Diaphragm and Filter			●
Illumination system	Plan-concave Mirror	●	●	●
	Attachable Illuminator	○	○	○
Filter	Blue Filter			●
	Green Filter			●

Note: ● Standard Outfit, ○ Optional