

BS-2030MH Multi-Head Microscope



BS-2030MH4A



BS-2030MH4B



BS-2030MH10

Introduction

BS-2030MH Series Multi-Head Microscopes are equipped with multi-head for more persons observing at the same time. The optical system is high quality and reliable.

Feature

1. Sharp image display with excellent optical system.
2. Small space occupancy with integral stand design, low environment requirement with anti-mould technology.
3. Users friendly and comfortable operation with ergonomic design.

Application

BS-2030MH Series Multi-Head Microscopes are mainly used in medical teaching and biological teaching. They also can be used for biological analysis for more experts at the same time. They can be widely used in medical and sanitary establishments, laboratories, institutes, colleges and universities.

Specification

Item	Specification	BS-2030 MH4A	BS-2030 MH4B	BS-2030 MH10
Viewing Head	Sliding Binocular Head, inclined at 45°, Interpupillary Distance: 55-75mm, Anti-mould	2PCS	1PC	5PCS
	Monocular Head, inclined at 45°, Anti-mould		2PCS	
	Sliding Trinocular Head, inclined at 45°, Interpupillary Distance: 55-75mm, Anti-mould	○	○	○
Eyepiece	Wide Field Eyepiece WF10×/18mm	4PCS	4PCS	10PCS
	Wide Field Eyepiece WF16×/11mm			
Nosepiece	Quadruple Nosepiece	●	●	●
Objective	Achromatic Objective 4×, 10×, 40×, 100×	●	●	●
	Achromatic Objective 20×, 60×	○	○	○
	Semi- Plan Achromatic Objective 4×, 10×, 20×, 40×, 60×, 100×	○	○	○
	Plan Achromatic Objective 4×, 10×, 20×, 40×, 60×, 100×	○	○	○
Condenser	Abbe Condenser NA1.25	●	●	●
Focusing	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 28mm	●	●	●
Stage	Double Layers Mechanical Stage 140×140mm/ 75×50mm	●	●	●
Illumination	1W S-LED illumination, Brightness adjustable	●	●	●
	Halogen Lamp 6V/ 20W	○	○	
	Halogen Lamp 24V/ 100W			○
Pointer	Green LED Pointer, Brightness Adjustable			●

Note: ●Standard Outfit, ○Optional