

BS-3003 Long Working Distance Stereo Microscope



Introduction

BS-3003 Long working distance stereo microscope is a smart and cost-effective stereo microscopes. A variety of eyepieces and objectives are available for different requirements. The long working distance feature provides more choice for different samples and works. It is widely used in electric factories, schools laboratories, sculpture, families and so on.

Feature

1. Sharp stereo erect image over a wide field of view.
2. 45° Inclined binocular head with locked-in eyepiece. The left ocular-tube with diopter adjustment $\pm 5dp$.
3. V-shape inserted objective with locked-in structure, easy to replace.
4. Come with high brightness 3W LED lamp.
5. A wide range of accessories for selection.

Application

BS-3003 stereo microscope is widely applied in school education, circuit board inspection and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts.

Specification

Item	Specification	BS-3003
Viewing Head	Binocular Viewing Head, 45° Inclined, Interpupillary Distance 55-75mm, diopter adjustment $\pm 5dp$ on left ocular-tube, Tube diameter 23.2mm	●
Eyepiece	WF5×/20mm	○
	WF10×/20mm	●

	WF15x/15mm	○
	WF20x/10mm	○
Objective	0.5x, Working distance 250mm	○
	1x, Working distance 230mm	●
	2x, Working distance 80mm	○
Illumination	Incident 3W LED illumination	●
Power supply	110V/60HZ or 220V/50HZ, 110V-240V Wide input voltage is optional	●
Package	1pcs/carton, 41cm*20cm*43cm, Gross weight 7.5kg, Net weight 6kg	●

Note: ● Standard Outfit, ○ Optional