

## BSL-150A Series Cold Light Illumination



Cold Light Power Supplier



### Introduction

BSL-150A Series Cold Light Illumination is designed as an auxiliary lighting device for microscope in order to get better observation results.

### Feature

1. More stable power supply with CE standard electrical parts.
2. With delay startup design to protect eyes and electrical parts from sudden brightness.
3. Enhancing the security with skillful power off switch.
4. Cooling the illumination better with advanced wind loop design.
5. Reliable with stable structure.
6. Longer working life and lower noisy.

### Specification

Item	Specification	BSL-150A-1	BSL-150A-2	BSL-150A-O
Power Supplier	Switch Power 110V/ 220V	●	●	●
	T-H Lamp with Reflector 20V/ 150W			
	Brightness Adjustable			
	Multi-coated IR Cutoff Filter for Heat Insulation			
	Fan Cooling			
	Size: 260×120×150mm			
	Weight 3.3 kg			
Single Light Guide	Rigid Flexional Tube, Length 500mm	●		
Dual Light Guide	Rigid Flexional Tube, Length 500mm		●	
Ring Light Guide	Flexible Tube, Length 1000mm, Adapter Ring Size Φ50mm/ Φ60mm			●
Kohler Illumination Attachment	Attached on Single or Dual Light Guide to adjust and obtain Even Illumination	○	○	

Note: ● Standard Outfit, ○ Optional