

Blood Smear Microscope Slides



FEATURE

- * The unique hydrophilic surface is the key condition of successful blood smear.
- * High quality glass material, special cleaning process, precise cutting process and good specification consistency which greatly improving the work efficiency.
- * The ColorCoat surface can be used for traditional label identification, 2B pencil and designated marker, also ideal for use with inkjet, barcode and QR code slide marking machine.
- * Automatic Blood Smear Microscope Slides Specially developed for sample maker equipment in the market to make automatic blood smear slides, such as Sysmex fully-automated slide-maker Sp 1000i and BECKMAN COULTER LH755 fully-automated slide-maker etc.
- * Manual Blood Smear Microscope Slides is suitable for manual blood smear preparations, also ideal raw material for liquid-based cytology, e.g., TCT & LCT slide preparation.

SPECIFICATION

Item No.	Description	Dimension	Edges	Corner	Packaging	Color
RM7310A	Automatic	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	white, orange, green, pink, blue and yellow
RM7320A	Manual	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	white, orange, green, pink, blue and yellow

When ordering this model, please indicate the color.