

# Circular Cover Glass



## FEATURE

- \* Excellent optical properties, stable molecular structure, flat surface and highly consistent size.
- \* Recommended for manual workflow in histology, cytology, urinalysis and microbiology.

## SPECIFICATION

Item No.	Dimension	Thickness	Packaging	Glass
RM6208	Ø8	0.13-0.16mm	100pcs/box	Normal Glass
RM6208A	Ø8	0.13-0.16mm	100pcs/box	Super White Glass
RM6208B	Ø8	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6210	Ø10	0.13-0.16mm	100pcs/box	Normal Glass
RM6210A	Ø10	0.13-0.16mm	100pcs/box	Super White Glass
RM6210B	Ø10	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6211	Ø11	0.13-0.16mm	100pcs/box	Normal Glass
RM6211A	Ø11	0.13-0.16mm	100pcs/box	Super White Glass
RM6211B	Ø11	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6212	Ø12	0.13-0.16mm	100pcs/box	Normal Glass
RM6212A	Ø12	0.13-0.16mm	100pcs/box	Super White Glass

RM6212B	Ø12	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6213	Ø13	0.13-0.16mm	100pcs/box	Normal Glass
RM6213A	Ø13	0.13-0.16mm	100pcs/box	Super White Glass
RM6213B	Ø13	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6214	Ø14	0.13-0.16mm	100pcs/box	Normal Glass
RM6214A	Ø14	0.13-0.16mm	100pcs/box	Super White Glass
RM6214B	Ø14	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6215	Ø15	0.13-0.16mm	100pcs/box	Normal Glass
RM6215A	Ø15	0.13-0.16mm	100pcs/box	Super White Glass
RM6215B	Ø15	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6216	Ø16	0.13-0.16mm	100pcs/box	Normal Glass
RM6216A	Ø16	0.13-0.16mm	100pcs/box	Super White Glass
RM6216B	Ø16	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6218	Ø18	0.13-0.16mm	100pcs/box	Normal Glass
RM6218A	Ø18	0.13-0.16mm	100pcs/box	Super White Glass
RM6218B	Ø18	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6220	Ø20	0.13-0.16mm	100pcs/box	Normal Glass
RM6220A	Ø20	0.13-0.16mm	100pcs/box	Super White Glass
RM6220B	Ø20	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6222	Ø22	0.13-0.16mm	100pcs/box	Normal Glass
RM6222A	Ø22	0.13-0.16mm	100pcs/box	Super White Glass
RM6222B	Ø22	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6224	Ø24	0.13-0.16mm	100pcs/box	Normal Glass
RM6224A	Ø24	0.13-0.16mm	100pcs/box	Super White Glass
RM6224B	Ø24	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6225	Ø25	0.13-0.16mm	100pcs/box	Normal Glass
RM6225A	Ø25	0.13-0.16mm	100pcs/box	Super White Glass
RM6225B	Ø25	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6230	Ø30	0.13-0.16mm	100pcs/box	Normal Glass
RM6230A	Ø30	0.13-0.16mm	100pcs/box	Super White Glass
RM6230B	Ø30	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6237	Ø37	0.13-0.16mm	100pcs/box	Normal Glass
RM6237A	Ø37	0.13-0.16mm	100pcs/box	Super White Glass
RM6237B	Ø37	0.13-0.16mm	100pcs/box	Boro 3.3 Glass

## OPTIONAL

Different thickness optional to meet the personalized needs of different customers

Thickness: 0.09mm-0.13mm 0.16mm-0.19mm 0.13mm-0.16mm 0.19mm-0.23mm