


ColorCoat Microscope Slides



FEATURE

- *Pre-cleaned, ready for use.
- *Ground edges and 45° corner design which is greatly reduces the risk of scratching during the operation.
- * ColorCoat Slides come with a light opaque coating in six standard colors: white, orange, green, pink, blue and yellow, resistant to common chemicals and routine stains that are utilized in laboratory
- * One-sided paint, it will not change color in routine H&E staining.
- * Suitable for marking with inkjet and thermal transfer printers and permanent markers

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Category	Color
RM7109	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade	white, orange, green, pink, blue and yellow
RM7109A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade	white, orange, green, pink, blue and yellow

When ordering this model, please indicate the color.

OPTIONAL

Other options to meet the personalized needs of different customers.

Dimension	Thickness	Edges	Corner	Packaging	Category
25x75 mm	1-1.2mm	Ground Edges	45°	50pcs/box	Standard Grade
25.4x76.2mm (1"x3")		Cut Edges	90°	72pcs/box	Super Grade
26x76mm		Beveled Edges			