

## D type Diagnostic Microscope Slides

### FEATURE

\* Different wells are coated with PTFE according to the needs of customers. Due to the excellent hydrophobic property of PTFE coating, it can ensure that there is no cross contamination between the wells, which can detect multiple samples on a diagnostic slide, save the amount of reagent used, and improve the detection efficiency.

\* It is suitable for all kinds of immunofluorescence experiments, especially for the immunofluorescence disease detection kit, which provides an excellent solution for microscope slide.

### SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Marking surface	Additional coating	Wells
<b>RM7410D</b>	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	white	No coating	Multiple optional

When ordering this model, please indicate the aperture.

1 well, $\Phi$ 6mm 	1 well, $\Phi$ 8mm 	2 wells, $\Phi$ 8mm, with numbers 
2 wells, $\Phi$ 11mm 	3 wells, $\Phi$ 11mm 	3 wells, $\Phi$ 14mm 
4 wells, $\Phi$ 6mm, with numbers 	4 wells, $\Phi$ 11mm, with numbers 	5 wells, $\Phi$ 8mm, with numbers 
6 wells, $\Phi$ 5mm, with numbers 	8 wells, $\Phi$ 6mm, with numbers 	10 wells, $\Phi$ 6mm, with numbers 
12 wells, $\Phi$ 5mm, with numbers 	14 wells, $\Phi$ 5mm, with numbers 	18 wells, $\Phi$ 5mm, with numbers 