

## DSLR Camera Eyepiece Adapter

These 2 adapters are used to connect DSLR camera to the microscope eyepiece tube or trinocular tube of 23.2mm. If the eyepiece tube diameter is 30mm or 30.5mm, you can plug the 23.2 adapter into the 30mm or 30.5mm connecting ring and then plug into the eyepiece tube.



BDPL-1(NIKON)



BDPL-2(CANON)

	Camera Mount	Magnification	Connecting Dia.	Application
BDPL-1(NIKON) adapter	Nikon	2×	23.2mm	Used to connect Nikon DSLR camera to eyepiece / trinocular tube
BDPL-2(CANON) adapter	Canon	2×	23.2mm	Used to connect Canon DSLR camera to eyepiece / trinocular tube