


Frosted Microscope Slides



FEATURE

- *Pre-cleaned, ready for use.
- *Ground edges and 45° corner design which is greatly reduces the risk of scratching during the operation.
- * Frosted area is even and delicate, and resistant to common chemicals and routine stains that are utilized in laboratory
- * Meet most of the experimental requirements, such as histopathology, cytology and hematology, etc.

SPECIFICATION

Item No.	Frosted Side	Dimension	Edges	Corner	Packaging	Category
RM7105	Single Frosted	25x75, 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade
RM7105A	Single Frosted	25x75, 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade
RM7107	Double Frosted	25x75, 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade
RM7107A	Double Frosted	25x75, 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade

OPTIONAL

Other options to meet the personalized needs of different customers.

Frosted Side	Dimension	Thickness	Edges	Corner	Packaging	Category
Single Frosted	25x75 mm	1-1.2mm	Ground Edges	45°	50pcs/box	Standard Grade
Double Frosted	25.4x76.2mm (1"x3") 26x76mm		Cut Edges Beveled Edges	90°	72pcs/box	Super Grade