


Hydrophilic Adhesion Microscope Slides



FEATURE

- *Treated with several coating technologies, which make the slides have strong adhesion and hydrophilic surface.
- * Optimized for use with Roche Ventana IHC automated stainer.
- * Recommended for manual IHC staining, automatic IHC staining with the Dako, Leica and Roche Ventana IHC automated stainer.
- * Ideal for use in H&E staining for routine and frozen sections like fat section, brain section and bone section where require stronger adhesion.
- * Suitable for marking with inkjet and thermal printers and permanent markers.
- * Six standard colors: white, orange, green, pink, blue and yellow, which is convenient for users to distinguish different types of samples and alleviate the visual fatigue in work.

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Category	Color
RM7204	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade	white, orange, green, pink, blue and yellow
RM7204A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade	white, orange, green, pink, blue and yellow

When ordering this model, please indicate the color

OPTIONAL

Other options to meet the personalized needs of different customers

Dimension	Thickness	Edges	Corner	Packaging	Category
25x75 mm	1-1.2mm	Ground Edges	45°	50pcs/box	Standard Grade
25.4x76.2mm (1"x3")		Cut Edges	90°	72pcs/box	Super Grade
26x76mm		Beveled Edges			