

L type Diagnostic Microscope Slides

FEATURE

* Different wells are coated with PTFE according to the needs of customers. Due to the excellent hydrophobic property of PTFE coating, it can ensure that there is no cross contamination between the wells, which can detect multiple samples on a diagnostic slide, save the amount of reagent used, and improve the detection efficiency.

* Ideal for liquid-based slide preparation.

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Marking surface	Additional coating	Wells
RM7420L	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	white	No coating	Multiple optional

When ordering this model, please indicate the aperture.

1 well, $\Phi 13\text{mm}$



2 wells, $\Phi 13\text{mm}$



8 wells, $\Phi 8\text{mm}$, with numbers



OPTIONAL

Additional coating: no coating, polysine coating

Customized design and hydrophobic printing color are available.