


Liquid-Based Cytology Microscope Slides



FEATURE

- * Supplied for liquid-based cytology, e.g., TCT & LCT slide preparation.
- * The hydrophilic surface makes the cells spread more evenly on the surface of the slide, without a large number of cells stacking and overlapping. The cells are clearly visible and easy to observe and identify.
- * Suitable for marking with inkjet and thermal printers and permanent markers.
- * Six standard colors: white, orange, green, pink, blue and yellow, which is convenient for users to distinguish different types of samples and alleviate the visual fatigue in work

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Category	Color
RM7205	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade	white, orange, green, pink, blue and yellow
RM7205A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade	white, orange, green, pink, blue and yellow

When ordering this model, please indicate the color

OPTIONAL

Other options to meet the personalized needs of different customers

Dimension	Thickness	Edges	Corner	Packaging	Category
25x75 mm	1-1.2mm	Ground Edges	45°	50pcs/box	Standard Grade
25.4x76.2mm (1"x3")		Cut Edges	90°	72pcs/box	Super Grade
26x76mm		Beveled Edges			