


Microscope Slides with Cavity



FEATURE

- * Designed for examining living micro-organisms, such as bacteria and yeast in hanging drops.
- * Ground edges and 45° corner design which is greatly reduces the risk of scratching during the operation.

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Cavity No.	Cavity Dia.	Cavity Depth
RM7103A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Single	15-18mm	0.3-0.5mm
RM7104A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Double	15-18mm	0.3-0.5mm

OPTIONAL

According to customized customer requirements, thicker glass slides can be used to provide deeper cavity and ensure the consistency.