


Plain Microscope Slides



FEATURE

- *Pre-cleaned, ready for use.
- *Ground edges and 45° corner design which is greatly reduces the risk of scratching during the operation.
- *Recommended for routine H&E stains and microscopy in laboratory, can also be used as teaching experiments.

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Category
RM7101	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade
RM7101A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade

OPTIONAL

Other options to meet the personalized needs of different customers.

Dimension	Thickness	Edges	Corner	Packaging	Category
25x75 mm	1-1.2mm	Ground Edges	45°	50pcs/box	Standard Grade
25.4x76.2mm (1"x3")		Cut Edges	90°	72pcs/box	Super Grade
26x76mm		Beveled Edges			