


Polysine Adhesion Microscope Slides



FEATURE

- * Polysine Slide is pre-coated with Polysine which improve the adhesion of tissues to the slide.
- * Recommended for routine H&E stains, IHC, ISH, frozen sections and cell culture.
- * Suitable for marking with inkjet and thermal transfer printers and permanent markers.
- * Six standard colors: white, orange, green, pink, blue and yellow, which is convenient for users to distinguish different types of samples and alleviate the visual fatigue in work.

SPECIFICATION

Item No.	Dimension	Edges	Corner	Packaging	Category	Color
RM7202	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Standard Grade	white, orange, green, pink, blue and yellow
RM7202A	25x75mm 1-1.2mm Thick	Ground Edges	45°	50pcs/box	Super Grade	white, orange, green, pink, blue and yellow

When ordering this model, please indicate the color

OPTIONAL

Other options to meet the personalized needs of different customers

Dimension	Thickness	Edges	Corner	Packaging	Category
25x75 mm	1-1.2mm	Ground Edges	45°	50pcs/box	Standard Grade
25.4x76.2mm (1"x3")		Cut Edges	90°	72pcs/box	Super Grade
26x76mm		Beveled Edges			