

Square and Rectangular Cover Glass



FEATURE

- * Excellent optical properties, stable molecular structure, flat surface and highly consistent size.
- * It can be widely used in routine laboratory and pathology laboratory for histopathology, cytology, urinalysis, microbiology, etc.

SPECIFICATION

Item No.	Dimension	Thickness	Packaging	Glass
RM6118	18X18mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6118A	18X18mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6118B	18X18mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6120	20X20mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6120A	20X20mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6120B	20X20mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6122	22X22mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6122A	22X22mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6122B	22X22mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6124	24X24mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6124A	24X24mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6124B	24X24mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6132	24X32mm	0.13-0.16mm	100pcs/box	Normal Glass

RM6132A	24X32mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6132B	24X32mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6140	24X40mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6140A	24X40mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6140B	24X40mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6150	24X50mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6150A	24X50mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6150B	24X50mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass
RM6160	24X60mm	0.13-0.16mm	100pcs/box	Normal Glass
RM6160A	24X60mm	0.13-0.16mm	100pcs/box	Super White Glass
RM6160B	24X60mm	0.13-0.16mm	100pcs/box	Boro 3.3 Glass

OPTIONAL

Different thickness optional to meet the personalized needs of different customers

Thickness: 0.09mm-0.13mm 0.16mm-0.19mm 0.13mm-0.16mm 0.19mm-0.23mm