

Stand and other accessories



<p><b>BSZ-F1 Stand</b> Column height: 400mm, <math>\Phi</math>38mm Column for focus Arm: <math>\Phi</math>30mm</p>	<p><b>BSZ-F2 Stand</b> Column height: 454mm, <math>\Phi</math>29.5mm Base size: 396*276mm Column for focus Arm: <math>\Phi</math>30mm</p>	<p><b>BSZ-F3 Stand</b> Column height: 490mm, <math>\Phi</math>38mm Base size: 253*253mm, Cross bar length: 446mm, Column for focus Arm: <math>\Phi</math>30mm</p>
--	---	---



<p><b>BSZ-F4 Stand</b> Column for focus Arm: <math>\Phi</math>30mm</p>	<p><b>BSZ-F9 Stand</b> Column height: 280mm Glass plate: <math>\Phi</math>100mm Microscope mount: <math>\Phi</math>76mm</p>	<p><b>BSZ-F10 Stand</b> Column height: 280mm Glass plate: <math>\Phi</math>140mm Microscope mount: <math>\Phi</math>76mm</p>
--	---	--



<p><b>BSZ-F11 Stand</b> Column height: 280mm Glass plate: <math>\Phi</math>100mm Microscope mount: <math>\Phi</math>76mm</p>	<p><b>BSZ-F12 Stand</b> Halogen incident and transmitted illumination Microscope mount: <math>\Phi</math>76mm</p>	<p><b>BSZ-F13 Stand</b> Clamp Size: 103*77*115mm Max. Table thickness: 52mm Length of Horizontal Arm: 507mm Microscope mount: <math>\Phi</math>76/<math>\Phi</math>40/<math>\Phi</math>45mm</p>
--	---	---



<p><b>BSZ-F14 Stand</b> Base Size: 382*260*24mm Column height: 408mm, <math>\Phi</math>32mm Microscope mount: <math>\Phi</math>76/<math>\Phi</math>40/<math>\Phi</math>45mm</p>	<p><b>BSZ-F15 Stand</b> Base Size: 382*260*24mm Column height: 408mm, <math>\Phi</math>32mm Cross bar Length: 408mm Microscope mount: <math>\Phi</math>76/<math>\Phi</math>40/<math>\Phi</math>45mm</p>	<p><b>BSZ-F16 Stand</b> Base Size: 318*308*16mm Column height: 326mm Focusing Range: 160mm Microscope mount: <math>\Phi</math>76/<math>\Phi</math>40/<math>\Phi</math>45mm</p>
---	---	--



<p align="center"><b>BSZ-F17 Stand</b></p> <p>Base Size: 318*308*16mm          Column height: 326mm          Focusing Range: 160mm          Microscope Arm: <math>\Phi 76/\Phi 40/\Phi 45</math>mm</p>	<p align="center"><b>BSZ-F18 Stand</b></p> <p>Base Size: 318*308*16mm          Column height: 500mm          Focusing Range: 160mm          Microscope Arm: <math>\Phi 76/\Phi 40/\Phi 45</math>mm</p>	<p align="center"><b>BSZ-F19 Stand(For Camera)</b></p> <p>Base Size: 318*308*16mm          Column height: 326mm          Focusing Range: 160mm</p>
--	--	--



<p align="center"><b>BMS-302 X-Y Stage</b></p> <p>Stage Size: 180×155×26mm,          Moving Rang: X-75mm, Y-55mm,          Accuracy: 0.1mm</p>	<p align="center"><b>BMS-304 X-Y Stage</b></p> <p>Stage size: 254×203×25mm,          Moving Range: 100mm×100mm,          Glass Stage: <math>\Phi 140</math>mm</p>	<p align="center"><b>BMS-311 Measuring Stage</b></p> <p>Stage Size: 160×123×21mm,          Moving Rang: X-75mm, Y-50mm, Accuracy:          0.1mm</p>
--	---	--